

2025 Centrifuge Bundles

Offers valid until Dec. 31st, 2025

Z206-A



Z206-A Compact Bundle

- 6 x 50ml rotor



6,000 rpm / 4,427 xg

Z206-BND: **\$2,800**

Z216-MK



Z216-M Microcentrifuge Bundles

- 44 x 01.5/2.0 ml
- COMBI-rotor, 24 x 1.5/2.0ml and 4 x PCR
- Over 21,000 xg! (15,000 rpm)



44 x 1.5/2.0ml

24 x 1.5/2.0ml and 4-PCR strips

Ambient
Z216-BND44: **\$3,850**
(44-place)

Z216-CMB: **\$4,400**
(COMBI-rotor)

Refrigerated
Z216-K-BND44: **\$6,895**
(44-place)

Z216-K-CMB: **\$7,400**
(COMBI-rotor)



Z287-A



Z287-A MicroBundle

- 2 x microplate and 18 x 1.5/2.0ml



3,800 rpm / 3,717 xg

14,000 rpm / 15,994 xg

Z287-BND-MB: **\$5,300**

Z307



Z307 Deepwell Plate Bundle

- 2 x Deepwell Plate rotor



4,500 rpm / 2,716 xg

Z307-BND-MP: **\$6,000**

Z327-K



Z327 Series Bundles

- Tissue culture or blood tube options
- up to 8 x 50ml and 16 x 15ml
- up to 28 x 8-10ml and 32 x 1.8-7ml



Tissue Culture Pkg.



Blood Tube Pkg.

5,000 rpm / 4,136 xg

Ambient
Z327-BND: **\$9,200**
(Tissue Culture)

Z327-BND-BD: **\$9,200**
(Blood tube)

Refrigerated
Z327-K-BND: **\$12,200**
(Tissue Culture)

Z327-K-BND-BD: **\$12,200**
(Blood tube)

2025 Centrifuge Bundles

Offers valid until Dec. 31st, 2025

Z366-K



Z366 Series Bundles

- Ambient or Refrigerated
- up to 16 x 50ml and 40 x 15ml



4,500 rpm / 3,735 xg

Ambient
Z366-BND: **\$9,900**

Refrigerated
Z366-K-BND: **\$12,650**

Z446-K



Z446 Series Bundles

- Tissue culture or Microplate options
- up to 28 x 50ml and 68 x 15ml
- up to 20 std. plates and 4 x Deepwell



Ambient:
4,000 rpm / 3,434 xg

Refrigerated:
4,500 rpm / 4,346 xg

Tissue Culture Pkg. Microplate Pkg.

Ambient
Z446-BND: **\$12,200**
(Tissue Culture)
Z446-MPBND: **\$12,000**
(Microplate)

Refrigerated
Z446-K-BND: **\$15,700**
(Tissue Culture)
Z446-K-MPBND: **\$15,600**
(Microplate)



Z496-K



Z496-K Floor Standing Bundle

- Swing out rotor: up to 40 x 50ml and 96 x 15ml
- Refrigerated



4,200 rpm / 3,964 xg

Only available in 230V

Z496-K-BND: \$28,600

Z36-HK



Z36-HK Super Speed Bundles

- Max speed up to 65,395
- 3 rotor options available
- Refrigerated



6 x 250ml rotor
10,000 rpm / 15,650 xg

6 x 50ml rotor
21,000 rpm / 41,410 xg

12 x 1.5/2.0 ml rotor
30,000 rpm / 65,395 xg

Z36-BND250: \$17,700
(6 x 250 ml)

Z36-BND50: \$16,500
(6 x 50 ml)

Z36-BND20: \$15,500
(12 x 1.5/2.0 ml)

Z600-S3



SIEVA III Filtration Bundle

- Includes perforated basket *filters sold separately
- One-of-a-kind benchtop filtration centrifuge
- Option with temp. equalization port



Perforated rotor
10,000 rpm / 7,825 xg

*Available filter pore sizes: 7µ, 20µ, 30µ, 50µ

Z600-S3-BND: \$22,100

Z600-S3-TBND: \$23,100
(w/ temp. equalization)