

AURELIA® GLOVES
A SUPERMAX HEALTHCARE BRAND

PERFORMANCE

YOU CAN MEASURE



- ✓ Chemo-tested by the FDA
- ✓ Precision grip for lab work
- ✓ Consistent quality, reliable supply
- ✓ Multiple options for different lab needs

FDA 510k Approved Examination Gloves



- Sold in 100 and 300 count options
- 3.2mil finger thickness
- Cobalt Blue
- Finger textured



- Sold in 200 count boxes
- 3.5mil finger thickness
- Pink
- Finger textured



- Sold in 200 count boxes
- 3.5mil finger thickness
- Teal Green
- Finger textured



- Sold in 100 count boxes
- 4mil finger thickness
- White
- Finger textured



- Sold in 200 count boxes
- 4mil finger thickness
- Bluple
- Finger textured



- Sold in 100 count boxes
- 4.5mil finger thickness
- Sky Blue
- Finger Textured



- Sold in 100, 90, and 80 count boxes depending on size
- 4.5mil finger thickness
- Patriot Blue
- Finger Textured



Tested and Approved for Fentanyl & Chemo drug resistance per ASTM D6978-05

Authorized Aurelia Distributor:

midsci.com
custserv@midsci.com

800.227.9997
f: 636.225.9998

MIDSCITM
Laboratory Equipment & Supplies

AURELIA[®] GLOVES
www.aureliagloves.com