

# DeepDip™ Pipette Tips - 1500µL

*Fit for All Popular 1000, 1250, and 1300L Pipettes*

- ▶ Fits 1000µL universal single & multi-channel pipettors
- ▶ Reaches the bottom of conical centrifuge tubes
- ▶ Replaces 1000µL, 1250µL, and 1300µL tips

MTC Bio's DeepDip™ pipette tips are made for working in long, narrow openings such as 15mL and 50mL centrifuge tubes. They protect your pipettor from liquid and vapor contamination, allowing you to safely aspirate liquid samples without submerging the pipette shaft and protecting it from contact with the sidewalls.

These tips are 130mm long and fit on all popular 1000µL and 1300µL universal pipetting instruments (single and multi-channel) including Gilson®, Eppendorf® and many others.

They are e-beam sterilized and packaged in racks of 96 tips.



**The only tip that goes the distance!**

Cat. No.	Description	Qty
<b>DeepDip™ Extended Length Pipette tips</b>		
P4330	Pipette Tip, 1500µL Capacity, 130mm Length, Non-Filtered, Sterile, 5 Racks of 96 Tips	480/cs
P4331	Pipette Tip, 1500µL Capacity, 130mm Length, Filtered, Sterile, 5 Racks of 96 Tips	480/cs