

Exact Volume Transfer Pipettes


Sterile Packages, an MTC Exclusive

- ◉ Dual-bulb design ensures consistent volume delivery
- ◉ Manufactured from chemically resistant LDPE
- ◉ Minimal liquid retention
- ◉ Available in 10 capacities
- ◉ **STERILE** Resealable bags



Designed for consistency and ease of use, MTC Bio Exact Volume Transfer Pipettes deliver repeatable liquid handling without the complexity of calibrated instruments. Engineered for clean release and dependable performance, they are ideal for routine lab workflows, molecular applications, and kit-based use where simplicity and reliability matter most.

Overdraw Bulb
Captures excess liquid to deliver consistent and replicable volumes



The upper bulb facilitates pipetting with a familiar pinch and release, while the lower overdraw bulb captures excess liquid. This extra sample is trapped in the overdraw bulb during release, ensuring that only the desired volume is dispensed in pipetting.

E9500-20
20µL Capacity
53mm Length



E9500-50
50µL Capacity
72mm Length



E9500-80
80µL Capacity
96mm Length



E9500-100
100µL Capacity
86mm Length



E9500-150
150µL Capacity
100mm Length



E9500-200
200µL Capacity
100mm Length



E9500-250
250µL Capacity
113mm Length



E9500-300
300µL Capacity
113mm Length



E9500-400
400µL Capacity
115mm Length



E9500-500
500µL Capacity
120mm Length



Cat. No.	Description	QTY
E9500-20	Pipette, Dual Bulb, 20µL, 53mm, Sterile, 10 Resealable Bags of 500	500
E9500-50	Pipette, Dual Bulb, 50µL, 72mm, Sterile, 5 Resealable Bags of 500	500
E9500-80	Pipette, Dual Bulb, 80µL, 96mm, Sterile, 5 Resealable Bags of 500	500
E9500-100	Pipette, Dual Bulb, 100µL, 86mm, Sterile, 5 Resealable Bags of 500	500
E9500-150	Pipette, Dual Bulb, 150µL, 100mm, Sterile, 5 Resealable Bags of 500	500
E9500-200	Pipette, Dual Bulb, 200µL, 100mm, Sterile, 5 Resealable Bags of 500	500
E9500-250	Pipette, Dual Bulb, 250µL, 113mm, Sterile, 5 Resealable Bags of 500	500
E9500-300	Pipette, Dual Bulb, 300µL, 113mm, Sterile, 2 Resealable Bags of 500	500
E9500-400	Pipette, Dual Bulb, 400µL, 115mm, Sterile, 2 Resealable Bags of 500	500
E9500-500	Pipette, Dual Bulb, 500µL, 120mm, Sterile, 2 Resealable Bags of 500	500