

DESIGNED TO ENHANCE WORKER PERFORMANCE.



PF-95BL

POWDER-FREE ECOTEK BIODEGRADABLE
NITRILE EXAMINATION GLOVES WITH DRY CORE*

**KEEP HANDS COOLER AND MORE COMFORTABLE
WITH PREMIUM EXAM QUALITY NITRILE GLOVES,
FEATURING UNIQUE DRY CORE MOISTURE
MANAGEMENT TECHNOLOGY.**

- Dry Core® Moisture Management Technology pulls sweat away from the skin
- EcoTek® Biodegradability Technology accelerates landfill breakdown without any performance loss*
- Silicone-free for critical environments

PF

PF-95BL

WEIGHT	LENGTH	THICKNESS	FEATURES	COLOR
 8.5g (large size)	 9.5"	 6.1 mil	Dry Core® EcoTek®	 BLUE


ECOTEK



ANSI Level 3
Abrasion
Protection

Fully textured
surface for wet
and dry grip

Tear-resistant
beaded cuff for
easy donning





PF-95BL

POWDER-FREE ECOTEK BIODEGRADABLE
NITRILE EXAMINATION GLOVES WITH DRY CORE*

PRODUCT DESCRIPTION

SPECIFICATIONS

PRODUCT		PowerForm®	
GRADE		Powder-Free Exam	
MATERIAL		Nitrile	
THICKNESS	FINGER	6.9mil	0.18mm
	PALM	6.1mil	0.16mm
	CUFF	5.1mil	0.13mm
COLOR	OUTER	Blue	
	INNER	White	
SURFACE	OUTER	Fully Textured	
	INNER	Dry Core®	
LENGTH		9.5"	242mm
WEIGHT		8.5g	
FIT TYPE		Ambidextrous	
CUFF TYPE		Beaded	
WATERTIGHT AQL		1.5	
ORIGIN COUNTRY		China	

PHYSICAL PROPERTIES

	PRE-AGING	POST-AGING
TENSILE STRENGTH:	20 MPa	19 MPa
ELONGATION:	560%	550%

INTEGRATED TECHNOLOGIES

DRY CORE Moisture Management Technology: Absorbent lining wicks away moisture for added comfort in extended use.

ECOTEK Biodegradability Technology: 92.6% biodegradable in 2.5 years. ASTM D5526-12 test results do not predict future biodegradation.*



* Specifications based on size Large glove

STANDARDS & CERTIFICATIONS

EN ISO 21420:2020
EN 455-2000-2015
ASTM D6319-19
ASTM F1671-13
ASTM F739-12
Reach EU 1907/2006
EU 2016/425

ANSI/ISEA 105-2016
CAT III Complex Design
FDA 21 CFR 177.2600



Please visit swsafety.com/symbols for information on standards used on this document.

ORDERING

ITEM NO.	SIZE	PALM WIDTH	PACKAGING
N106551	XS 5½ – 6	3.0" 76mm	50 gloves per box 10 boxes per case
N106552	S 6½ – 7	3.4" 87mm	50 gloves per box 10 boxes per case
N106553	M 7½ – 8	3.8" 97mm	50 gloves per box 10 boxes per case
N106554	L 8½ – 9	4.3" 109mm	50 gloves per box 10 boxes per case
N106555	XL 9½ – 10	4.5" 114mm	50 gloves per box 10 boxes per case
N106556	2XL 10½ – 11	4.9" 125mm	50 gloves per box 10 boxes per case



APPLICATIONS

INDUSTRIAL

Agriculture
Hydroponics
Assembly

MEDICAL

Chemical Handling
Clinical
Laboratory
Liquid Handling
Pharma Production
Tattooing and Cosmetics

Keep out of sunlight. Store in a cool, dry place
Keep away from sources of ozone or ignition



R_N10655X1121