

White (W) is ideal for gPCR. White or colored plates do not effect regular PCR runs.

Skirt Type			non-Skirted Plates		semi-Skirted					full-Skirted			Corbett
Profile / Well Volume			0.1 mL	0.2 mL	0.1 mL	0.1 mL	0.2 mL	0.2 mL	0.2 mL	0.1 mL	384	384	Corbett 0.1 mL
Thermal Cycler Brand	Cycler Model	Block Format	PR-PCR1196 PR-PCR1196-W	PR-PCR1296 PR-PCR1296-W PR-PCR1296-RB PR-PCR1296-CUT	PR-PCR2196F	PR-PCR2196LC PR-PCR2196LC-W	PR-PCR2296	PR-PCR2296R	PR-PCR2296S PR-PCR2296S-W	PR-PCR3196 PR-PCR3196-W	PR-PCR33841 PR-PCR33841-W	PR-PCR33842	PR-PCRCorbett
ABI Applied Biosystems	2700, 2720, 7000, 7300, 7500, 9700, MiniAmp, MiniAmp Plus, Proflex, Quant Studio 3, Quant Studio 5, Quant Studio 6, Quant Studio 7, QuantStudio 12k, SimpliAmp, Veriti, Veriti Dx, ViiA7, 7900 HT	96 x 0.2mL					x	x	x				
	7500, 9800, Quant Studio 3, Quant Studio 5, Quant Studio 6, Quant Studio 7, QuantStudio 12k, Veriti, Veriti Dx, ViiA7, 7900 HT, StepOne Plus	96 x 0.1mL Fast/Low Profile			x								
	9700, Automated Thermal Cycler, Proflex, Quant Studio 5, Quant Studio 6, Quant Studio 7, QuantStudio 12k, Veriti, Veriti Dx, ViiA7, 7900 HT	384									x		
	Automated Thermal Cycler	96 x 0.2mL						x					
Agilent	Gradient Cycler, RoboCycler	96 x 0.2mL	x	x					x	x			
	Mx3000P, Mx3005P, Mx4000	96 x 0.2mL		x					x				
	SureCycler 8800	96 x 0.2mL	x	x									
	AriaMx	96 x 0.1mL Fast/Low Profile	x							x			
Analytik Jena / Biometra	Gradient Cycler, SureCycler 8800	384									x		
	FlexCycler2, qTower, SpeedCycler2, T1, TAdvanced, TGradient, TOne, TProfessional, TRobot	96 x 0.2mL	x	x					x	x			
	T1, Tadvanced, Tgradient, Trobot	384									x		
	Axygen	MaxyGene II	96 x 0.2mL		x				x				
Bioer / Bulldog	GeneTouch, GeneTouch Plus, LifeECO, LifePro, LifeTouch	96 x 0.2mL	x	x					x	x			
	GeneTouch	384									x		
Bio-Rad	C1000/C1000 Touch, Chromo4, DNA Engine, MJ Mini, PTC-200, S1000	96 x 0.2mL	x	x					x	x			
	iCycler, iCycler IQ, IQS, MyIQ/MyIQ2, T100	96 x 0.2mL		x					x				
	CFX Connect, CFX Touch, Opticon	96 x 0.1mL Fast/Low Profile	x							x			
	C1000/C1000 Touch, CFX Touch, DNA Engine, PTC-200, S1000	384									x		
Bruker / Hain Lifesciences	FluoroCycler XT	96 x 0.1mL Fast/Low Profile				x							
	GTQ-Cycler 96	96 x 0.2mL	x	x					x	x			
Cole Parmer / Techne	TC-S12, TC-PLUS	96 x 0.2mL	x	x					x	x			
	Flexigene, Genius, Prime, Prime Elite, TC-4000, TC-412, TC-5000,	96 x 0.2mL		x					x				
	PrimeQ, Quantica,	96 x 0.1mL Fast/Low Profile	x							x			
	Flexigene, Prime, Prime Elite, TC-4000, TC-412, TC-5000, TC-PLUS	384									x		
DLab	TC1000-G	96 x 0.2mL	x	x					x	x			
Eastwin	ETC811, ETC821	96 x 0.2mL		x					x				
Eppendorf	Mastercycler ep Gradient, Mastercycler ep realplex, Mastercycler Gradient, Mastercycler Nexus, Mastercycler Pro, Mastercycler X50	96 x 0.2mL	x	x					x	x			
	Mastercycler ep Gradient, Mastercycler Gradient, Mastercycler Pro, Mastercycler X50	384									x		
Esco	Aeris, Swift MaxPro	96 x 0.2mL	x	x					x	x			
	Aeris, Swift MaxPro	384									x		
Genestorm / Agilegene	GS1, GS2, GS4	96 x 0.2mL	x	x					x	x			
	GS1, GS2, GS4	384									x		
Kyratec	SuperCycler Trinity	96 x 0.2mL	x	x					x	x			
Labnet	MultiGene/ OptiMax	96 x 0.2mL		x					x				
MWG Biotech	Primus, The Q LifeCycler	96 x 0.2mL	x	x					x	x			
	Primus	384									x		
PegLab / VWR	peqSTAR 2X, peqSTAR X, Uno	96 x 0.2mL	x	x					x	x			
	peqSTAR X, Uno	384									x		
Qiagen	Rotor-Gene Q	4 x 0.1mL											x
Roche®	LightCycler 480, LightCycler 96	96 x 0.1mL				x						x	
SensoQuest	Labcycler	96 x 0.2mL	x	x					x	x			
	Labcycler	384									x		
Takara	TP 3000	96 x 0.2mL	x	x					x				
ThermoFisher	Arktik, Hybaid MBS, Hybaid Px2/PxE	96 x 0.2mL	x	x					x	x			
	Arktik, Hybaid MBS, Hybaid Px2/PxE	384									x		