

- Volume adjustment with plunger button or thumbwheel
- Volume easily visible with black and red numbers
- Reduced ejection forces
- Fully autoclavable for sterility
- Adjustable ejector length
- Pipettors are individually tested and supplied with certificate of quality
- Shafts and tip cone with higher strength and chemical resistance
- Accurate and Precise
- AlphaPette™ Pipettor have a universal shaft that will accept most pipet tips, unlike Rainin LTS Pipettors



AlphaPette™ Pipettors are designed for exceptional comfort and accuracy. Combining ease of use: lightweight; well balanced; comfortable fit in the right or left hand; and accuracy, the AlphaPette™ is the perfect addition to your lab.

The adjustable thumbwheel is slightly recessed to prevent unintentional change and is protected from the effects of thermal expansion by the insulated body of the pipette. For convenience, volume can also be adjusted by turning the push-button.

AlphaPette™ Pipettors have a universal shaft that will accept most pipette tips. The stainless-steel tip ejector is adjustable to accommodate various styles of tips and removable for pipetting into narrow tubes. Ease of use meets exceptional performance with the AlphaPette™.

## PR1MA™

ITEM #	DESCRIPTION	QUANTITY
A-2	0.1 - 2 µL	1/ea
A-10	0.5 - 10 µL	1/ea
A-20	2 - 20 µL	1/ea
A-100	10 - 100 µL	1/ea
A-200	20 - 200 µL	1/ea
A-1000	100 - 1000 µL	1/ea
A-5000	1000 - 5000 µL	1/ea