

- Superior accuracy and reproducibility
- Lightweight and balanced for comfort
- Lower plunger forces reduce strain
- Curved ejector bar lessens ejection forces
- Universal shaft will accept most pipette tips
- Tip ejector is adjustable and removable
- Continuously adjustable Vol selection with thumbwheel or plunger
- Fully autoclavable
- 3-Year Warranty

These precision pipettes are designed for exceptional comfort and accuracy. Lightweight and well-balanced, they fit comfortably in either the right or left hand with the large digital display visible to the user. Chemically resistant, maintenance-free seals and finely polished pistons provide unsurpassed accuracy sample after sample.

The continuously adjustable, digital micrometer is slightly recessed to prevent unintentional change and is protected from the effects of thermal expansion by the insulated body of the pipette. For convenience, Vol can also be adjusted by turning the push button. Ergonomically positioned, the tip ejector reduces fatigue during prolonged use.

The AlphaPette™ Plus Pipettes have a universal shaft that will accept most pipette tips. The tip ejector is adjustable to accommodate various styles of tips and is removable for pipetting into narrow tubes. Made from stainless steel, it will not rust, corrode, or break. The AlphaPette™ Plus Pipettes are fully autoclavable.



PR1MA™

ITEM #	DESCRIPTION	QUANTITY
AP4508-10	8-Channel, 1 to 10 µL	1/ea
AP4508-50	8-Channel, 5 to 50 µL	1/ea
AP4508-200	8-Channel, 20 to 200 µL	1/ea
AP4508-300	8-Channel, 50 to 300 µL	1/ea
AP4512-10	12-Channel, 1 to 10 µL	1/ea
AP4512-50	12-Channel, 5 to 50 µL	1/ea
AP4512-200	12-Channel, 20 to 200 µL	1/ea
AP4512-300	12-Channel, 50 to 300 µL	1/ea