

- Keep PCR Tubes, Strip Tubes, and PCR Plates (96 & 384) cold or hot on the bench
- Non-Skid rubber feet provide a stable platform for loading samples
- Gets cold in as few as 25 minutes
- Samples stay cold on your bench top for up to an hour
- Small footprint minimizes bench top and freezer space used
- Precisely drilled wells allow for maximum heat transfer between the block and samples
- Small footprint minimizes benchtop and freezer space used
- Precisely drilled wells allow for maximum heat transfer between the block and samples



ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
RACK-CRYO	15 x Cryo Tubes, for Cryo Tubes, Silver	1/pk
RACK-96	96 x 0.2 mL, for PCR, Silver	1/pk
RACK-96-RB	96 x 0.2 mL, for PCR, Assorted	1/pk
RACK-384	384-Plates, for PCR, Silver	1/pk
RACK-15	15 x 1/5/2.0 mL, for Microcentrifuge Tubes, Silver	1/pk
RACK-5	24 x 0.5 mL, for Microcentrifuge Tubes, Silver	1/pk