



- Ceramic work surface, 7.1 x 7.1 in.
- Stir speed from 60 to 1,500 rpm
- Temperature control up to 380°C
- Three models: heat & stir, heat only or stir only



PR1MA™ hotplates, stirrers, and hotplate-stirrers feature an exceptionally durable, chemical-resistant white ceramic work surface. Their space-efficient design (8 x 9 in. footprint), makes them ideal for use on crowded benchtops and inside of bio-hoods.

Advanced microprocessor controls with convenient turn knobs allow quick, precise adjustment and maintenance of speed and temperature. Safety indicator LED's on the front panel indicate when the heating and/or stirring functions have been activated. With a square 7.1 inch work surface, all three models

are compatible with a wide variety of popular sizes of borosilicate glass beakers, flasks, bottles, and other vessels. The hotplate/magnetic stirrer model (MS-1500) comes complete with a support rod for mounting thermometers and temperature probes.

SPECIFICATIONS

SPEED RANGE	60 to 1,500 rpm
TEMP. RANGE	80° C to 380° C
DIMENSIONS	8.5 x 10 x 4.5 in. 20 x 23 x 11.5 cm
WEIGHT	9 lbs 4kg
ELECTRICAL	120 V, 60Hz
MAX STIRRING CAPACITY	3.5 Liters
WARRANTY	2-Year Warranty

ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
HP-380	Hotplate 7.1 x 7.1 in., 115 V	1/ea
MS-1500	Magnetic Stirrer 7.1 x 7.1 in., 115 V	1/ea
SH-1500	Hotplate & Magnetic Stirrer 7.1 x 7.1 in., 115 V	1/ea

