



- Digital temperature control
- Personal sized, 13 x 14.5 in. footprint (20 L)
- Accepts bottles and flasks up to 2 L
- Internal outlet, for installation of shaker, rockers, etc.

With digital temperature control, the PR1MA™ Digital Refrigerated Incubator provides a convenient “set and walk away” operation, eliminating the need for external thermometers and repetitive “fine-tuning” of an analog control knob. Simply choose the desired temperature and the incubator immediately begins to heat up while accurately monitoring the chamber temperature and conveniently displaying the temperature in real-time on the large LED control panel.

Despite a modest footprint of 14.5 x 13 in., the incubator features a large internal chamber, capable of storing flasks and bottles up to 2 Liters. In addition, the incubator includes two adjustable shelves and an internal outlet for powering small shakers, rockers, and other instrumentation. A mini nutating rocker (H3D1020) supplied with both flat and dimpled mats is available for the agitation of samples during incubation.

## SPECIFICATIONS

|                            |   |
|----------------------------|---|
| <b>TEMP. RANGE</b>         | 5°C above ambient to 60°C                     |
| <b>TEMP. ACCURACY</b>      | 0.5°C (at 37°C)                               |
| <b>TEMP. UNIFORMITY</b>    | ± 1.5°C (at 37°C)                             |
| <b>CAPACITY</b>            | 20 L   0.75 ft <sup>3</sup>                   |
| <b>EXTERIOR DIMENSIONS</b> | 13.2 x 14.5 x 18.7 in.<br>33.5 x 37 x 47.5 cm |
| <b>WEIGHT</b>              | 15 lbs   6.5 kg                               |
| <b>ELECTRICAL</b>          | 100 to 240V, 50-60 Hz                         |
| <b>WARRANTY</b>            | 2-Year Warranty                               |

## ORDERING INFORMATION

| ITEM #   | DESCRIPTION                                     | QUANTITY |
|----------|---|----------|
| IN20-H   | Heat Only (Ambient +5 to 60°C), Incl. 2 Shelves | 1/ea     |
| IN20-HC  | Mini Digital Incubator - Heat and Cool          | 1/ea     |
| H2200-SH | Extra Shelf, Plastic                            | 1/ea     |