

# PR1MA™ Mini Dry Baths

# PR1MA™



- Available in heat only or heat & cool versions
- Exchangeable blocks for 0.2 to 50 mL tubes
- Simple touchpad, digital control
- Link two consecutive programs (cooling model only)
- Custom blocks available upon request

With its compact construction, broad temperature range, and multiple block options, the PR1MA™ Mini is the perfect tool for laboratory incubations. The simple touchpad controls and a digital display provide for "set and walk away" temperature selection and unrivaled accuracy.

While perfect for bench use, the optional 12V adapter expands the utility of the PR1MA™ Mini to field use. At less than 4.3 inches wide, PR1MA™ Mini is truly the first personal block incubator.

## SPECIFICATIONS

MODEL	MB300	MB100
TEMP RANGE	RT -25°C to 100°C	RT +5°C to 100°C
DIMENSIONS	4.3 x 6.4 x 5.5 in. 11 x 16.2 x 14 cm	4 x 6 x 4.5 in. 10.5 x 15 x 11.5 cm
WEIGHT	2 lbs   1 kg	
ELECTRICAL	Universal Voltage, 100 to 240 V, 50-60 Hz, 50 W	
WARRANTY	2-Year Warranty	



## ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
MB100	Mini Digital Dry Bath	1/ea
MB300	Mini Digital Dry Bath, Heat and Cool	1/ea

Blocks sold separately

**MIDSCI™**  
Laboratory Equipment & Supplies



## PR1MA™ Mini Digital Dry Bath with Heated Lid (MB100-HL)

- Heated lid prevents evaporation & condensation
- Compact, fits in the palm of your hand
- Independent temperature control of lid & block

The PR1MA™ Mini Digital Dry Bath with Heated Lid is the first and only mini heat block featuring a heated lid for protecting samples against unwanted condensation and evaporation.

Both the block and lid temperature can be selected independently in increments of 0.1°C, making it the ideal choice for the most temperature-sensitive applications, including immunoassays, LAMP assays, enzyme reactions & other temperature hyper-sensitive procedures.

### SPECIFICATIONS

<b>TEMP RANGE</b>	Ambient +5 to 100°C
<b>DIMENSIONS</b>	4 x 6 x 4.5 in. 10.5 x 15 x 11.5 cm
<b>WEIGHT</b>	2 lbs   1 kg
<b>ELECTRICAL</b>	Universal Voltage, 100 to 240V, 50-60Hz, 50W
<b>WARRANTY</b>	2-Year Warranty

### ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
MB100-HL	Mini Digital Dry Bath, Heated Lid	1/ea