

- Digital speed control up to 1,500 rpm
- Compact, 3 mm orbit for thorough mixing in microplates
- Includes platform for up to 4 microplates
- Cold room and incubator safe
- Optional adapters for mixing microtubes



The PR1MA™ Orbital Microplate Shaker is designed with a 3 mm mixing orbit, optimized for thorough mixing of microplates and/or PCR plates. Up to four plates (standard or deep well) can be placed on the included platform and are instantly secured without the use of springs, clamps, or tools.

## SPECIFICATIONS

<b>SPEED</b>	Variable, 200 to 1,500 rpm
<b>MAX CAPACITY</b>	Up to 4 microplates
<b>DIMENSIONS</b>	13.1 x 11.74 x 4.7 in. 33 x 30 x 12 cm
<b>WEIGHT</b>	12 lbs   5.5 kg
<b>ELECTRICAL</b>	Universal Voltage, 100-240 V, 50-60 Hz 250 W
<b>WARRANTY</b>	2-Year Warranty

## ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
BT1500	13.1 x 11.7 x 4.7 in. Shaker, Holds Up to 4 Microplates	1/ea
BT1500-A05	20 x 0.5 mL Optional Adapter	1/ea
BT1500-A1520	20 x 1.5/2.0 mL Optional Adapter	1/ea