


PR1MA-Shell Rigid PCR Plates

- Rigid polycarbonate frame to reduce warping commonly associated with polypropylene PCR plates.
- Ultra-thin polypropylene wells provide optimal thermal transfer and higher reaction efficiencies.
- White well options enhance the fluorescence signal during qPCR.
- Broad thermal cycler compatibility to fit most commercially-available instruments.
- Raised rim around each well for improved sealing.
- High contrast printed alphanumeric.
- ANSI/SBS footprint and well spacing make them perfectly suitable for liquid handling instruments.
- Certified DNase, RNase, and pyrogen-free.

Online PCR
Consumables
Selector



Item #	# of Wells	Volume	Description	Frame Color	Well Color	Packaging	List Price
PR-PCR3742-C	96	0.1 mL	Full-Skirt	Clear*	Clear	10 plates/pk	\$53.00
PR-PCR3742-CW	96	0.1 mL	Full-Skirt	Clear	White	10 plates/pk	\$61.00
PR-PCR3750-C	96	0.2 mL	Semi-Skirt	Clear*	Clear	10 plates/pk	\$53.00
PR-PCR3735-C	384	30 µL	Full-Skirt, Single Notch	Clear*	Clear	10 plates/pk	\$75.00
PR-PCR3735-CW	384	30 µL	Full-Skirt, Single Notch	Clear	White	10 plates/pk	\$79.00
PR-PCR3730-C	384	30 µL	Dual-Notch	Clear	Clear	10 plates/pk	\$75.00

* Available in individual colors (blue, green, red, violet, yellow) or in an assortment of the five colors.

PR1MA Polypropylene PCR Plates

- Broad thermal cycler compatibility to fit most commercially-available instruments.
- Raised rim around each well for improved sealing.
- Optically clear with printed black alphanumeric.
- Ultra-thin polypropylene wells provide optimal thermal transfer and higher reaction efficiencies.
- Certified DNase, RNase, and pyrogen-free.




Item #	# of Wells	Volume	Description	Color	Notch	Packaging	List Price
PR-PCR1296	96	0.2 mL	No-Skirt	Clear	A12	10 plates/pk	\$27.00
PR-PCR1196	96	0.1 mL	No-Skirt, Low-Profile	Clear	H12	10 plates/pk	\$53.00
PR-PCR2296S	96	0.2 mL	Semi-Skirt	Clear	A12	10 plates/pk	\$27.00
PR-PCR2296	96	0.2 mL	Semi-Skirt, ABI®-Type	Clear	A12	10 plates/pk	\$27.00
PR-PCR2296R	96	0.2 mL	Semi-Skirt, ABI®-Type, Raised Rim	Clear	A12	10 plates/pk	\$27.00
PR-PCR2196F	96	0.1 mL	Semi-Skirt, Low Profile, Raised Sides, ABI Fast®-Type	Clear	A1	10 plates/pk	\$27.00
PR-PCR2196LC-W	96	0.1 mL	Semi-Skirt, Low Profile, LightCycler® Compatible	White	H12	10 plates/pk	\$32.00
PR-PCR3196	96	0.1 mL	Full-Skirt, Low-Profile	Clear	A12	10 plates/pk	\$27.00
PR-PCR3196-W	96	0.1 mL	Full-Skirt, Low-Profile, Bio-Rad® Compatible	White	A12	10 plates/pk	\$32.00
PR-PCR33841	384	30 µL	Full-Skirt, ABI®-Type	Clear	A24	10 plates/pk	\$39.00
PR-PCR33842	384	30 µL	Full-Skirt, Roche®-Type	Clear	A24, P24	10 plates/pk	\$32.00



PCR Strip Tubes



Combo Pack
PR-PCR28CF



Rigid Strips feature robust, **3-point interwell linkages** making strips more sturdy and decreasing the chance of spillage.



Break-Away Strips feature single interwell linkages, making strips easy to split into smaller sections.



Optically-clear caps ensure that the fluorescence signal is passed through uninhibited.

Item #	Type	Volume	Description	Packaging	List Price
PR-PCR28ACF	Rigid	0.2 mL	8-Strip Tubes with Individually Attached Hinge Caps	120 strips/pk	\$73.00
PR-PCR18ACF	Rigid	0.1 mL	8-Strip Tubes with Individually Attached Hinge Caps, Low-Profile	120 strips/pk	\$75.00
PR-PCR18ACF-W	Rigid	0.1 mL	8-Strip Tubes with Individually Attached Hinge Caps, Low-Profile, White	120 strips/pk	\$78.00
PR-PCR28ACD	Break-Away	0.2 mL	8-Strip Tubes with Individually Attached Dome Caps	120 strips/pk	\$73.00
PR-PCR28ACFF	Break-Away	0.2 mL	8-Strip Tubes with Individually Attached Flat Caps	120 strips/pk	\$73.00
PR-PCR28CD	Combo Pack	0.2 mL	8-Strip Tubes and Dome Caps	125 strips/pk	\$81.00
PR-PCR28CF	Combo Pack	0.2 mL	8-Strip Tubes and Flat Caps	125 strips/pk	\$81.00
PR-PCR18CF	Combo Pack	0.1 mL	8-Strip Tubes and Flat Caps, Low-Profile	125 strips/pk	\$81.00

Individual PCR Tubes

Item #	Volume	Description	Packaging	List Price
PR-PCR2D	0.2 mL	Dome Cap PCR Tube	1,000/pk	\$48.00
PR-PCR2F	0.2 mL	Optically-Clear, Flat Cap PCR Tube	1,000/pk	\$41.00
PR-PCR2FF	0.2 mL	Frosted, Flat Cap PCR Tube	1,000/pk	\$41.00



Online PCR Consumables Selector

Your MIDSCI Sales Representative:

Name: _____

Email: _____

Phone: _____