

- Accepts most popular pipettes
- Stable, sturdy steel base
- Rotates 360 degrees
- Six pipette capacity



Designed to accommodate a set of six single-channel pipettes and keep them easily in reach while also keeping your benchtop organized. The freely rotating circular rack is universally fitting for the vast majority of pipettes.

With a convenient handle knob at the top, the entire carousel stand can be easily and safely carried or repositioned, even with a full load of six pipettes. Keep your pipettes safe with the wide heavyweight base that contains a steel stability plate that eliminates tipping. Compact, stable, and versatile, this carousel stand is a welcome addition to any research laboratory.

ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
MID4280-CR	Universal Pipette Carousel, Holds 6 Pipettors	1/ea