

- Buy one that lasts; 2-year warranty
- 8-hour rechargeable battery life
- QuickStand™ retractable kickstand
- For pipettes from 1 mL to 100 mL

PR1MA™ Pipet Controllers are lightweight (only 6.3oz/170g) and designed to fit comfortably in either the right or the left hand. The natural feel of the ergonomic body design, the computer balanced weight distribution and the low-pressure fingertip control buttons... all combine to provide a more comfortable pipetting experience.

The attached QuickStand™, retractable kickstand, allows the controller to stand on the bench without removing or emptying the serological pipet. The stand also holds the pipet at an angle which prevents the flowing of liquid to the filter and nosepiece. When it's not being used, the QuickStand™ flips back into neutral position where it doesn't interfere with normal operation.

A long-life rechargeable lithium (environmentally friendly) battery keeps the controller running for up to 8 hours or 2500 pipettings on a single two-hour charge. PR1MA™ Pipette Controller's easy-fit silicone nosepiece accepts pipette sizes from 1 mL to 100 mL. With its light-weight advanced ergonomics, continuously variable electronic speed adjustment and extended lithium battery life, the PR1MA™ Pipette Controller is the logical choice for comfort, convenience and value.



## ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
MPIP-1	PR1MA™ Pipette Controller	1/ea
MID6950	PR1MA™ Reach Long Neck	1/ea