




Supreme accuracy with the low-retention feature!

- Barrier tips have extra working volume, so liquid doesn't touch the filter
- Great stackable racks
- 10 µL XL tips have 1 µL graduations
- Clear & readable graduations
- Sterile, RNase free and DNase free



Improve consistency and reproducibility with low-retention tips!

1 µL graduation on 10 µL XL tip

Graduations:

- 10 µL: 2 & 10 µL
- 10 µL XL: 1, 2, & 10 µL
- 20 µL: 10 µL
- 100 µL: 10, 50, & 100 µL
- 200 µL: 10, 50, & 100 µL
- 300 µL XL: 10, 50, 100, & 200 µL
- 1,000 µL: 100, 500, 1,000 µL
- 1,250 µL: 100, 500, & 1,000 µL

PR1MA™ Low-Retention Barrier Tips

ITEM #	DESCRIPTION	QUANTITY
PR-10RK-FL	10 µL, Sterile	96/rk, 10 rks/cs
PR-10XLRK-FL	10 µL XL, Sterile	96/rk, 10 rks/cs
PR-20RK-FL	20 µL, Sterile	96/rk, 10 rks/cs
PR-100RK-FL	100 µL, Sterile	96/rk, 10 rks/cs
PR-200RK-FL	200 µL, Sterile	96/rk, 10 rks/cs
PR-200XLRK-FL	200 µL XL 5 in. Long, Sterile	96/rk, 8 rks/cs
PR-300RK-FL	300 µL, Sterile	96/rk, 10 rks/cs
PR-1000RK-FL	1,000 µL, Sterile	96/rk, 8 rks/cs
PR-1250RK-FL	1,250 µL, Sterile	96/rk, 8 rks/cs

