



RMP2024



- Includes tube holders for common tube sizes, 11 x 1.5/0.2 mL (pack of 2), 5 x 10-15 mL (pack of 2), 5 x 5-7 mL (pack of 2), 6 x 50 mL (pack of 1)
- Plus model available with variable speed control and advanced shaking/ mixing modes: Mix (Oscillating), Pause, Rock



RFP2020

The new PR1MA™ Roto-Mini Series of rotating mixers provide a wide range of mixing capabilities. Suitable for an array of laboratory applications including resuspension, protein binding, nucleic acid purification, etc. the mixers include interchangeable tube holders for accommodating all common tube sizes from 0.5 mL to 50 mL. The motion can be adjusted from vertical (end over end) mixing, to horizontal (rolling or hybridization style) mixing.

SPECIFICATIONS

MODEL	RFP2020	RMP2024
SPEED	24 rpm	5 to 70 rpm
MIXING MODES	Mix, Pause, Rock	
DIMENSIONS	10.8 x 6 x 6.5 in 27.5 x 15.2 x 16.5 cm	
WEIGHT	3 lbs 1.4 kg	
ELECTRICAL	115 or 230 V, 50-60 Hz, 0.5A	Universal Voltage 100 to 240 V, 50-60 Hz, 0.5 A
WARRANTY	2-Year Warranty	

ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
RFP2020	10.8 x 6 x 6.5 in., 24 rpm, Fixed Speed	1/ea
RMP2024	10.8 x 6 x 6.5 in., 5 to 70 rpm, Variable Speed	1/ea

