

- All serological pipets are manufactured to the ASTM Standard Specification for Serological Pipets E-934
- Bi-directional (positive and negative) graduations are printed in bold easy-to-read ink
- Premium material barrier plug
- Color-coded stripes and packaging enable quick volume identification
- Printed indicator arrow for quickly locating negative graduations when using a pipet controller
- Pipets are made from Virgin Polystyrene
- Meets USP Class VI Certification Standards
- Animal Free
- Sterile, RNase/DNase Free, Non-pyrogenic: Validated to have less than 0.1 to 0.5 EU/mL



PR1MA™ Individually Wrapped

ITEM #	DESCRIPTION	QUANTITY
PR-SERO-1	1 mL, Individ. Wrapped, Sterile	800/cs
PR-SERO-2	2 mL, Individ. Wrapped, Sterile	600/cs
PR-SERO-5	5 mL, Individ. Wrapped, Sterile	200/cs
PR-SERO-10	10 mL, Individ. Wrapped, Sterile	200/cs
PR-SERO-25	25 mL, Individ. Wrapped, Sterile	200/cs
PR-SERO-50	50 mL, Individ. Wrapped, Sterile	90/cs

PR1MA™ Bulk Packed in Bags

ITEM #	DESCRIPTION	QUANTITY
PR-SERO-1SB	1 mL, Sterile	25 bulk packed/bag, 1000/cs
PR-SERO-2SB	2 mL, Sterile	25 bulk packed/bag, 1000/cs
PR-SERO-5SB	5 mL, Sterile	25 bulk packed/bag, 1000/cs
PR-SERO-10SB	10 mL, Sterile	25 bulk packed/bag, 1000/cs
PR-SERO-25SB	25 mL, Sterile	25 bulk packed/bag, 1000/cs

