

PR1MA™ Shell PCR Plates



Thermal Cycler Compatibility Guide

	96-Well Semi-Skirt 0.2mL	96-Well Full-Skirt 0.1mL		384-Well Full-Skirt A24 Notch		384-Well Full-Skirt A24, P24 Notches		
	Clear Wells	Clear Wells	White Wells	Clear Wells	White Wells	Clear Wells	White Wells	
APPLIED BIOSYSTEMS® (ABI)								
3100, 3130XL, 3500, 3500XL, 3730, 3730XL				●	●			
GeneAmp® 2700, 2720	●							
GeneAmp® 9600, 9700	●							
MiniAmp™, MiniAmp™ Plus	●							
ProFlex™	●							
Proflex™, Multiblock System				●	●			
QuantStudio™ 3, 5, 6, 7, 12 K Real-time PCR System (with 0.2mL blocks)	●			●	●			
SimpliAmp™	●							
Veriti™ 384-well				●	●			
Veriti™ 96-well (with 0.2mL block)	●							
ViiA7, 7900 HT FAST				●	●			
AGILENT®								
AriaMx qPCR System		●	●					
Mx3000P, Mx3005P, Mx4000	●							
Robocycler Gradient	●	●	●					
Surecycler 8800	●			●	●			
ANALYTIK JENA®								
Trobot 96, SpeedCycler2 SP, SPR								
Alpha SC								
Flexcycler2	●	●	●	●	●			
T1 Thermocycler, T Gradient, T One, T Advanced, T Professional series	●	●	●	●	●			
qTOWER3 (G, touch), Toptical	●			●	●			
BIOER® / BULLDOG®								
GeneExplorer™	●							
GeneTouch 384				●	●			
GeneTouch, GeneTouch Plus	●							
LifeECO™	●							
LifeTouch	●							
LifePro	●							



Every effort has been made to note registered trademarks. All trademarks, trade names, or logos mentioned or used are the property of their respective owners.



ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998

PR1MA™ Shell PCR Plates



Thermal Cycler Compatibility Guide

	96-Well Semi-Skirt 0.2mL	96-Well Full-Skirt 0.1mL		384-Well Full-Skirt A24 Notch		384-Well Full-Skirt A24, P24 Notches		
	Clear Wells	Clear Wells	White Wells	Clear Wells	White Wells	Clear Wells	White Wells	
GE / MEGABACE®								
GE Healthcare Amersham MegaBACE™ 4000				●	●			
GE Healthcare Amersham MegaBACE 500, MegaBACE™ 1000 mark 2		●	●					
ROCHE®								
LightCycler® 480 (384 well)						●	●	
TECHNE® / COLE-PARMER®								
Flexigene	●	●	●	●	●			
Techne® Genius, Genius Quad	●	●	●	●	●			
Techne PCRmax® Alpha cycler 1, 2, 4				●	●			
Techne® Prime, PrimeG, Prime Elite, Prime Elite Satellite	●			●	●			
Techne Quantica®		●	●	●	●			
Techne® TC412, TC512		●	●	●	●			
Techne® TC-PLUS	●	●	●					
Touchgene®, Touchgene® Gradient	●	●	●					
OTHER BRANDS / MODELS								
Apollo ATC-201	●							
BIONEER AllInOneCycler™	●							
DLab TC1000-G	●	●	●					
G-Storm® GS1, GS2, GS4	●	●	●					
HiMedia® Prima-96	●							
Labnet MultiGene™, OptiMax								
Idaho Technology LightScanner®		●	●	●	●			
peqLab peqSTAR 384X								
peqLab peqSTAR 96X	●							
Sensoquest Labcycler	●			●	●			
Takara Dice™ touch, Gradient	●							
ThermoFisher Hybaid MBS, Px2/PxE	●							
VWR® Collection UNO96	●	●	●					
VWR peqlab peqSTAR 96X		●	●					



Every effort has been made to note registered trademarks. All trademarks, trade names, or logos mentioned or used are the property of their respective owners.



ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998

PR1MA™ Shell PCR Plates



Thermal Cycler Compatibility Guide

	96-Well Semi-Skirt 0.2mL	96-Well Full-Skirt 0.1mL		384-Well Full-Skirt A24 Notch		384-Well Full-Skirt A24, P24 Notches		
	Clear Wells	Clear Wells	White Wells	Clear Wells	White Wells	Clear Wells	White Wells	
BIOMETRA®								
Biometra T Advanced				●	●			
Biometra T Advanced 96, 96 G, 96 S, 96 SG	●	●	●					
Biometra T One	●	●	●					
Biometra T Professional series	●	●	●					
Biometra Uno	●	●	●					
BIO-RAD® / MJ RESEARCH®								
C1000, C1000 Touch, S1000	●	●	●	●	●			
CFX Connect, CFX96, CFX96 Touch		●	●					
CFX384, CFX384 Touch				●	●			
CFX96 Touch Deep Well	●	●	●					
Chromo4™	●	●	●					
DNA Engine Opticon® and Opticon® 2		●	●					
DNA Engine, DNA Engine Tetrad®, Tetrad® 2		●	●	●	●			
DNA Engine Dyad®, Dyad® Disciple		●	●	●	●			
iCycler™	●							
iQ™, iQ™-4, iQ™-5, myiQ™ PCR System	●							
MJ Research PTC200, 220, 221, 225, 240	●			●	●			
MyCycler™	●							
PTC 200		●	●	●	●			
T100								
CORBETT® / QIAGEN®								
Corbett Research (Qiagen) Palm Cycler	●							
Corbett Research (Qiagen) Palm Cycler 384				●	●			
EPPENDORF®								
MasterCycler® ep 384 / Pro 384				●	●			
Mastercycler® ep realplex	●	●	●					
MasterCycler® ep, ep gradient, Pro, Pro S	●	●	●					
Nexus, Nexus gradient, Nexus SX1, Nexus GSX1	●	●	●					

Every effort has been made to note registered trademarks. All trademarks, trade names, or logos mentioned or used are the property of their respective owners.

PR1MA™

MIDSCI™
Laboratory Equipment & Supplies

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998