

PR1MA™ qMax Instrument Compatibility Table

REAL TIME INSTRUMENT BRAND	MODEL NUMBER/NAME	COMPATIBLE PR1MA QMAX™ REAGENTS	
Applied BioSystems®	7000, 7300, 7700 7900, 7900HT, 7900HT Fast StepOne™, StepOne™ Plus	qMax Green High Rox PR2000-H-100, PR2000-H-500, PR2000-H-1000, PR2000-H-S	qMax Probe High Rox PR2001-H-100, PR2001-H-500, PR2001-H-1000, PR2001-H-S
	7500, 7500 Fast ViiA7™ QuantStudio™ 12K Flex	qMax Green Low Rox PR2000-L-100, PR2000-L-500, PR2000-L-1000, PR2000-L-S	qMax Probe Low Rox PR2001-L-100, PR2001-L-500, PR2001-L-1000, PR2001-L-S
BioRad®	CFX96™ CFX384™	qMax Green No Rox PR2000-N-100, PR2000-N-500, PR2000-N-1000, PR2000-N-S	qMax Probe No Rox PR2001-N-100, PR2001-N-500, PR2001-N-1000, PR2001-N-S
	iCycler™, MyIQ™ IQ5™		
	Opticon™, Opticon™ 2 MiniOpticon™, Chromo4™	qMax Green No Rox PR2000-N-100, PR2000-N-500, PR2000-N-1000, PR2000-N-S	
Roche®	LightCycler™ 96 LightCycler™ 480 LightCycler™ Nano	qMax Green No Rox PR2000-N-100, PR2000-N-500, PR2000-N-1000, PR2000-N-S	qMax Probe No Rox PR2001-N-100, PR2001-N-500, PR2001-N-1000, PR2001-N-S
Qiagen® (Corbett)	Rotor-Gene™ Q Rotor-Gene™ 3000 Rotor-Gene™ 6000	qMax Green No Rox PR2000-N-100, PR2000-N-500, PR2000-N-1000, PR2000-N-S	qMax Probe No Rox PR2001-N-100, PR2001-N-500, PR2001-N-1000, PR2001-N-S
Eppendorf®	Mastercycler™ ep realplex Mastercycler™ ep realplex 2S	qMax Green No Rox PR2000-N-100, PR2000-N-500, PR2000-N-1000, PR2000-N-S	qMax Probe No Rox PR2001-N-100, PR2001-N-500, PR2001-N-1000, PR2001-N-S
Agilent® (Stratagene)	MX4000™ MX3000P™ MX3005P™ AriaMx	qMax Green Low Rox PR2000-L-100, PR2000-L-500, PR2000-L-1000, PR2000-L-Sx	qMax Probe Low Rox PR2001-L-100, PR2001-L-500, PR2001-L-1000, PR2001-L-S
Cepheid®	SmartCycler™	qMax Green No Rox PR2000-N-100, PR2000-N-500, PR2000-N-1000, PR2000-N-S	qMax Probe No Rox PR2001-N-100, PR2001-N-500, PR2001-N-1000, PR2001-N-S
Fluidigm®	BioMark™	qMax Green Low Rox PR2000-L-100, PR2000-L-500, PR2000-L-1000, PR2000-L-S	qMax Probe Low Rox PR2001-L-100, PR2001-L-500, PR2001-L-1000, PR2001-L-S
Illumina	Eco qPCR		qMax Probe No Rox PR2001-N-100, PR2001-N-500, PR2001-N-1000, PR2001-N-S

